# FIRE TESTING

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*EQUIPMENT FOR THE MAIN  
STANDARDS OF THE FIRE  
REACTION TESTING OF PLASTICS  
AND BUILDING MATERIAL*

ITA/ENG

[WWW.NOSELAB-ATS.COM](http://WWW.NOSELAB-ATS.COM)

***noselab* ats**  
ADVANCED TESTING SOLUTIONS

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**NOSELAB ATS**, *forty years of experience in the field of instrumentation for laboratory tests, in research and development of equipment to apply methods for the classification to reaction to fire of materials.*

*Over forty instruments and apparatuses have been developed and realized that allow to carry out the tests in a professional and repetitive way, guaranteeing to the companies and to the authorized certification laboratories, the validity of the process and the simplification of the test.*

*The guarantee is the continuous updating of the reference standards that deal with the analysis of the parameters of fire behavior of materials:*

**FLAMMABILITY**

**OXYGEN INDEX**

**SMOKE OPTICAL DENSITY**

**FLAME RESISTANCE**

**INCOMBUSTIBILITY**

**RADIANT EXPOSURE**

*These tests measure the parameters on a small scale, for other products where the conditions can not be adequately simulated an apparatus has been identified for the INTERMEDIATE ANALYSIS such as the CALORIMETRIC CONE.*

*The introduction of this apparatus is particularly important in the laboratories that have to deal with the major International Institutions all the indices are monitored by certified high precision instruments.*

STANDARD	MODELLO	EURO CLASS	BUILDINGS MATERIAL	AVIATION AUTOMOTIVE	POLYMERS TEXTILE	CABLES	SMOKE DENSITY
FAR 25 part 1	AA03			•			
DIN 75200 ISO 3795	AA04			•			
FAR 25 F part 1	AA05 AA07			•			
ISO 11925-2 DIN 4102	AA06	•	•				
IEC 60331-11-21 CEI 2036/1	AB01					•	
IEC 60332-1-1 CEI EN 50265-2-1	AB03		•			•	
IEC 60332-3-10	AB05		•			•	
ISO 9151 UNI EN 367	AC06				•		
ISO 6940 ISO 6941 ISO 15025	AC10				•		
UNI EN 9174	BA01	•	•				
ISO 9239 ASTM E 648	BA03	•	•				
ISO 5858 –2 IMO A 653	BA04			•			
ASTM E 662 ISO 5659-2 AITM 2.0007	CA01 CA03 CA06		•				•
CEI EN 61034-1	CA05		•			•	•
CEI 50267 –1 IEC 60754-1	DA01					•	
ISO 4589/2 LOI	EA04				•		
ISO 4589/3 LOI HT	EA05				•		
DIN 695 IEC 60695-2-10	EB03 EB05					•	
UL 94	ED01	•			•		
ISO 1182.2	FA01	•	•				
ISO 5660 ASTM E 1354	GA01	•	•	•	•		•

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*CABIN for fire reaction tests with horizontal sample, also available with timer to control flame application times.*

- Gas burner with BUNSEN diam. 9.5 mm (propane and methane nozzles)
- Gas regulation valve
- Extractable U-shaped sample holder, equipped with Ncr metal wire spacers and precipitated parts recovery pan
- 10091112 Model with Timer for checking the flame application time with automatic starter when inserting the sample, accuracy 0.1 sec

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*Allows you to test materials for construction and furniture. A flame is applied to the surface and to the lower edge of the specimen in a vertical position.*

*The parameters to be detected: Post-combustion time; Post-glow time; Detection of the damaged area; Emission of ignited drops.*

*The M5 burner for propane gas (LPG) equipped with microvalve to adjust the flame is mounted on a forward slide.*

*You can test various formats according to the dimensions required by the standards, using the different U-shaped holders.*

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
*The apparatus is useful for tests that rely on a flame with a controlled heat output of 750°C. Band burner for propane gas with three rows of holes, length 500 mm, Complete with Venturi mixer and fixing support, ignition device with safety devices From the console the regulation of the*

***The variants are available according to the types of cables and standard procedures.***

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*It allows you to determine the critical flux of radiant heat on the surface of materials cladding used on noncombustible substrates, especially for the production of walls. The radiation is applied by means of a gas powered ceramic panel. Duct for detecting hot smoke equipped with 5 thermocouples equipped with a system for recording the flame extension times along the surface of the vertically oriented sample subjected to radiant heat.*

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*CA01 Smoke Density Chamber in compliance with ASTM E 662 for determine the specific optical density of smoke generated by materials solids mounted in vertical position, with thicknesses up to 25 mm.*

*The use of a radiant energy source and a six burner flames reproduces irradiation conditions at the same time and combustion. While the smoke produced by the specimen accumulates in the sealed chamber a photometric system determines the optical transmittance variations. Model CA03 in accordance with the ISO 5659 standard adopts a furnace radiant heat (900W) to test in a horizontal position.*

*Photometric system with vertical light beam with parallel optics, silicon photodiode detector with spectral response similar to the S-4. Transmittance measurement range from 100 to 0.00001 in four scales, automatic zeroing.*

*Test data and specific optical density management are displayed on Personal Computers. DataLink Software handles them all the functions of opening and closing valves, calibration of the light beam, as well as to proceed with the test automatically. Reports of test with transmittance value diagrams, test archive.*

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*A gas sampling system can be adopted, FT-IR with gas chromatographic analysis using EGA modules for extraction, filtration and thermal stabilization of exhaust gases (up to 10 l / min). System developed according to EN 45545-2 , ISO 5659-2.*

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The test determines the flammability characteristics of insulating materials, sheaths and electrical cable fillers.

The oxygen index is the minimum volumetric oxygen percentage of a mixture of oxygen and nitrogen that maintains the combustion of the sample. The percentage of oxygen is determined by a paramagnetic transducer and the combustion takes place inside the column.

The detection of the oxygen percentage is given by a paramagnetic transducer (accuracy O<sub>2</sub> <0.1%).

Direct reading on the 7 "color display of % Oxygen, Flow, Air.

Automatic flow control through oxygen / nitrogen mixing massflow. Mixture production with addition of air for cheaper measures with high oxygen concentration. Automatic calibration (0 and 100%).

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EA05 with the same characteristics as the EA04 model, but to measure the minimum oxygen concentration required to support flame combustion in a mixture of oxygen and nitrogen at temperatures up to 400 ° C.

In addition to supplying the column power supply, the control unit manages the thermoregulator system for monitoring the column and sample temperature.

It has the control of the surface temperature of the specimen by means of a K-type thermocouple and thermoregulators for the temperature control of the preheating of the mixture - max 400 ° C and for the temperature control inside the glass column, PID type regulators.

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Alimentazione/Power supply: 230 V 50Hz

10095203 GLOW WIRE Automatico/ Automatic Glow wire

The appliance allows to simulate the overheating conditions, which materials such as electrical components or parts of insulation produced with thermoplastic materials could be subjected to malfunction or current overload. The electrically heated Glow Wire tip is pressed against the test sample at a defined temperature.

The system for the movement of the carriage allows to slide the specimen rack at a speed between 10 and 25 mm/sec.

Contact force against incandescent wire preloaded at  $0,95 \text{ N} \pm 0.1 \text{ N}$ .

All operations are set by color Touch Screen:

- wire temperature NiCr max  $960^\circ \text{C}$
- current in Ampere (True RMS)
- time of application of the incandescent wire contact time (TA) between test specimen and heating element ( $30 \pm 1$  seconds) and display of the ignition time (TI) and extinguishing (TE).

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
Apparatus with the same characteristics as the EB03 model but set up inside a steel test cabin of vol. about 0.5 cubic meters - internally finished in black to keep the level of lux controlled. The upper terminal is designed for connection to a smoke eva-cuation system, equipped with an internal light for control.

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
*Flammability test cabin according to UL94 (Underwriters Laboratories USA). The volume of about 1mc satisfies the required amount of oxygen to conduct the test even with extended testing times according to the following Test Methods:*

94HB ( ASTM D 635, ISO 1210) - Prova orizzontale/*Horizontal test*  
 94 V-0 V-1-V-2 (ASTM D 3801) - Prova verticale/*Vertical test*  
 94 5VA - 5VB (ASTM D 5048, ISO 9772) - Fiamma da 125 mm/ *flame 125 mm*  
 VTM-0-1-2 (ASTM D 4804, ISO 9773) - Prova verticale materiali non rigidi/ *Vertical test non-rigid materials*  
 HBF-HF 1-2 (ASTM D 4986) - Prova orizzontale materiali polimerici / *Horizontal test Polymeric material.*

Dimensioni/*Dimensions:*  
 1200x800x1700 h mm  
 Peso /*Weight:* 80 Kg

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*Ease of movement of the test accessories inside the booth, operation from the outside. Temperature measurement system and flame regulation (length 20 mm and 125 mm according to ASTM D 5207) with two thermocouples. Touch Screen control panel for flame application timing to detect postcombustion time and flame temperature reading.*

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
### Standard

ISO 1182

*Furnace in refractory tube heated by electrical resistors, thermally insulated with magnesium oxide; Diffuser cone insulated with fiberglass (0.04 W / mK); The positioning system is self-centering with manual movement and shock absorber.*

*The digital display for: temperature control up to 750 ° C, oven voltage regulation during the cycle (1.5kVA) alarm signal for exceeding the temperature.*

*Data acquisition system: Notebook and software with saving functions, data management, graphs, temperature and combustion time reports.*

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*Heat Release Rate is the key measure needed to assess the risk of fire of materials and products, quantifies the size of the fire, the rate of fire growth and consequently the release of toxic fumes and gases.*

*The release of the heat rate is used to classify the various products in terms of risk in the event of fire. The most accurate and effective device for measuring the release of heat is the Cone Calorimeter. The operating principle is based on the method of oxygen consumption, has found wide acceptance universally following the results usable for a range of fire models of furniture, cables, building materials, etc.*

*Several standardization bodies have published product standards for the use of the calorimeter to evaluate the performance of finished products:*

- ▣ Furniture (ASTM E 1474)
- ▣ Wall coating materials (ASTM E 1740)
- ▣ Mattress (ASTM F 1550)
- ▣ Electric cables (ASTM D6113)
- ▣ Rolling railway applications (BS 6873)
- ▣ Maritime applications (IMO)

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
### Functions:

- ▣ Heat release rate in kW / m<sup>2</sup>
- ▣ Real combustion heat in MJ / kg
- ▣ Mass loss speed
- ▣ Power-on time
- ▣ Density of fumes
- ▣ Soot sampling
- ▣ CO<sub>2</sub> and CO production
- ▣ Trigger time
- ▣ Other details related to combustion

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*The scintillator and the sample holder operated tire are automatically controlled by the software. Combustion is performed with a heat source of controlled radiation up to 100KW / m2 max. Temperature control is possible through the use of 3 K-type thermocouples and an electronic PID thermoregulator which also characterize the auto-tuning system to guarantee high stability.*

Gas analyzer:

- O<sub>2</sub>: paramagnetic
- 0-100% range and performance as per regulations
- <0.1% accuracy
- CO<sub>2</sub> and CO: sensors, non-dispersive infrared rays with performance in accordance with the standards
- 0-10% range for CO<sub>2</sub>
- 0-3.0% range for CO

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17 "TOUCH SCREEN VIDEO INTEGRATED FOR VISION AND THE MANAGEMENT OF TEST

*Combustion is carried out with a controlled radiation heat source up to 100KW/m2 max.*

*The load cell (0.05 g.) is mounted independently of the main panel for better stability and precision, thus avoiding any vibration caused by the exhaust fan. The measuring system is for high efficiency and equipped with protection by water cooling during the execution of the test.*

Functionality of the self-calibration system:

- Cell load
- GAS analysis
- Irradiance of the conical heater
- Smoke measurement system
- Management and display of status
- Control of the sensors and their relative stability
- Performance of all calculations required by standards (total and partial HRR, mass loss, etc.)
- All data can be acquired in CSV format to create excel files

***noselab* *ats***  
ADVANCED TESTING SOLUTIONS

NOSELAB ATS s.r.l.  
Via del Lavoro, 30  
20813 Bovisio Masciago (MB) - Italy  
Tel. +39 0362 367454 | +39 0362 450612  
info@noselab-ats.com  
www.noselab-ats.com

comimm.it

WWW.NOSELAB-ATS.COM